



## slimline3

The slimline alternating support system was designed with a view to balancing cost, features and effectiveness. Its twenty one high density lateral cells provide more touch points and therefore more patient comfort.

The advantage of not having static head cells and a CPR release means that the whole system can be rotated for positioning the control unit near the head of the bed.

Side handles for easy positioning

Quick release connectors

Anti kink air hoses

Heavy duty straps bond overlay to mattress

Use only with the consent and supervision of a suitably qualified health professional.  
 \*\* Mattress not included.

## product specification

Patient handling is greatly enhanced with the lower overlay profile and the heavy duty elastic straps and handles allow the overlay to be securely bonded to an underlying profiling mattress.\*\*

### Quick Facts

- Mattress overlay system
- Care level - Medium Risk
- Quick release zip on multi-stretch top cover
- Compact and quiet analogue control unit
- 21 high density cells alternating 1 in 2
- Transport feature
- 180 degree rotation for more choice in control unit positioning.

### control unit

|                           |                                       |
|---------------------------|---------------------------------------|
| control interface         | analogue                              |
| dimensions (WxHxD mm)     | 135x115x270mm                         |
| weight                    | 1kg                                   |
| air output                | 6.5 litres per minute                 |
| hose connection type      | quick release connectors              |
| timed static function     | <b>x</b>                              |
| alarms                    | visual                                |
| internal power down alarm | <b>x</b>                              |
| fuse rating (Amps)        | mains plug - 3A, control unit 1A (x1) |
| mode of operation         | continuous                            |
| electrical safety         | EN 60601.1:1996, UL2601-1 safe        |

### support surface

|                                |  |
|--------------------------------|--|
| support type                   | mattress overlay                         |
| dimensions (LxWxH millimetres) | 2000x880x115mm                           |
| alternation type               | 1 in 2                                   |
| alternation cycle time         | 12 minutes                               |
| base material                  | TPU coated nylon 420 denier woven fabric |
| cell material                  | TPU coated 210 denier nylon fabric       |
| cell in cell technology        | <b>x</b>                                 |
| side buffer                    | ✓ (optional)                             |
| number of cells                | 21                                       |
| static head cells              | <b>x</b>                                 |
| CPR function                   | <b>x</b>                                 |
| anti kink hosing               | ✓  |
| micro air loss feature         | <b>x</b>                                 |
| positioning handles            | ✓  |

### top cover

|                               |                               |
|-------------------------------|-------------------------------|
| top cover material            | PU coated 2 way stretch nylon |
| attachment type               | quick release zipper          |
| Infection control side skirts | <b>x</b>                      |
| seam joining                  | sewn                          |

### general

|                           |             |
|---------------------------|-------------|
| risk category             | medium risk |
| maximum user weight guide | 120kg       |
| carry bag included        | ✓           |